

What is new in imc LINK Version 1.2

Table of Contents

Changes and innovations in Version 1.2 Rev. 1

1.1 Windows 10 support	2
------------------------------	---

Changes and innovations in Version 1.2

2.1 Status Monitoring	2
2.2 Data transfer	2
2.3 Diskstart via FTP to internal device memory	2
2.4 imc FAMOS call	2
2.5 Logbook	2

Changes and innovations in Version 1.2 Rev. 1

1.1 Windows 10 support

imc LINK Version 1.2 Rev. 1 can be installed on Windows 10.

Changes and innovations in Version 1.2

2.1 Status Monitoring

Status Monitoring allows event-variables to be defined.

2.2 Data transfer

- Data can be retrieved from the device at adjustable time intervals (e.g. every 2 hours).
- Data can be retrieved from the device in response to changes event-variables.

2.3 Diskstart via FTP to internal device memory

Via FTP it is now possible to upload a configuration to the internal device memory.

2.4 imc FAMOS call

it is now possible to pass the device name (%5) and device identifier (%6) as parameters.

2.5 Logbook

Copying via context menu has been added.